



Lifespan Technologies
PO Box 3352
Salt Lake City, UT 84110
Tel: 801-971-1818 Fax: 801-464-6116
info@lifespantech.com

Product No: **K-1713**

**Low Molecular Weight Dextran
Sulfate Standards 6000-8000 MW**

K-1713: Low Molecular Weight Dextran Sulfate Standard Solutions

INTENDED USE: THIS PRODUCT IS FOR RESEARCH USE ONLY. NOT INTENDED FOR CLINICAL OR DIAGNOSTIC USE.

Description:

Low molecular weight (6000 – 8000 MW) dextran sulfate in tris buffered saline (TBS). Concentrations range from 0.0003 – 3 µg/ml and include a 0 µg/ml vial.

Storage and Stability:

Standards can be stored unopened and at 4°C for up to 6 months. If you are performing less than a full plate, the diluted standards may be stored tightly capped at 4° C and re-used with your next assay.

Standard Preparation:

The standards are supplied ready-to-use in Lifespan Technologies Dextran Sulfate ELISA kits. Mix well prior to pipetting.

Lifespan Technologies products are sold for research and development purposes only and are not for diagnostic use or to be incorporated into products for resale without written permission from Lifespan Technologies. Materials in this publication, as well as applications and methods and use, may be covered by one or more U.S. or foreign patents or patents pending. We welcome inquiries about licensing the use of our trademarks and technologies at MM@lifespantech.com.