



**Lifespan Technologies**  
PO Box 3352  
Salt Lake City, UT 84110  
Tel: 801-971-1818 Fax: 801-464-6116  
info@lifespantech.com

Product No: **K-1705**

**Unfractionated Heparin  
Standards**

## **K-1705: Unfractionated Heparin Standard Solutions**

**INTENDED USE: THIS PRODUCT IS FOR RESEARCH USE ONLY. NOT INTENDED FOR CLINICAL OR DIAGNOSTIC USE.**

### **Description:**

Unfractionated heparin in tris buffered saline (TBS).  
Concentrations range from 0.0003 – 30 µg/ml and include a 0 µg/ml vial.

### **Storage and Stability:**

Standards can be stored unopened and at 4°C for up to 6 months.  
If you are performing less than a full plate, the diluted standards may be stored tightly capped at 4° C and re-used with your next assay.

### **Standard Preparation:**

The standards are supplied ready-to-use in Lifespan Technologies Heparin ELISA kits.  
Mix well prior to pipetting.

Lifespan Technologies products are sold for research and development purposes only and are not for diagnostic use or to be incorporated into products for resale without written permission from Lifespan Technologies. Materials in this publication, as well as applications and methods and use, may be covered by one or more U.S. or foreign patents or patents pending. We welcome inquiries about licensing the use of our trademarks and technologies at [MM@lifespantech.com](mailto:MM@lifespantech.com).